



TERSO
SOLUTIONS, INC.

LARGE RFID CABINET TS056

- » REAL-TIME INVENTORY VISIBILITY
- » ENSURE PRODUCT INTEGRITY
- » ELIMINATE MANUAL PROCESSES

KNOW NOW.

Inventory Management Powered by RFID

BENEFITS OF THE TS056 - LARGE CABINET

- Optimize inventory levels
- Prevent stock-outs; monitor expiry of products
- Customize your product layout
- Provide regulatory compliance
- Monitor temperature remotely
- Protect critical inventory through secure access
- Tailor with your branding
- Integrate with existing inventory platforms



| COMPONENT | SPECIFICATION |
|-------------------------|--|
| External Dimensions | 2032 H x 914 W x 660 D (mm) |
| Internal Dimensions | 1651 H x 356 W x 610 D (mm) for each side |
| External Space Required | 2032 H x 1041 W x 1270 D (mm) Place 76mm from back wall Allow 534mm for door clearance |
| Volume | 717L |
| Weight | 120kg |
| Access Control | LF badge |

| COMMUNICATIONS | SPECIFICATION |
|---------------------|--|
| RAIN RFID Frequency | 864 to 868 MHz |
| Communications | Ethernet, optional cellular or wifi |

| POWER | SPECIFICATION |
|-------------------|------------------|
| Electrical Power | 230VAC, 50Hz, 1A |
| Power Consumption | 84W max |
| Battery Backup | Up to 1hr |