



TERSO
SOLUTIONS, INC.



LARGE RFID REFRIGERATOR TS074

- » REAL-TIME INVENTORY VISIBILITY
- » PROACTIVE TEMPERATURE ALERTING
- » ELIMINATE MANUAL PROCESSES

KNOW NOW.

Inventory Management Powered by RFID

BENEFITS OF THE TS074 - LARGE REFRIGERATOR

- Maintain temperature range of 2°C to 10°C
- Prevent stock-outs; monitor expiry of products
- Customize product layout
- Protect critical inventory through secure access
- Provide regulatory compliance
- Tailor with your branding
- Integrate with external temperature monitoring systems



COMPONENT	SPECIFICATION
External Dimensions	79”H x 30”W x 33”D
Internal Dimensions	58.3”H x 24.8”W x 24.3”D
External Space Required	87”H x 30”W x 66”D Place 3” from back wall Allow 21” for door clearance
Volume	20.2ft ³
Weight	500lbs
Cooling Method	Forced air circulation
Access Control	LF badge

COMMUNICATIONS	SPECIFICATION
RAIN RFID Frequency	902 to 928 MHz
Communications	Ethernet, optional cellular or wifi

POWER	SPECIFICATION
Electrical Power	115VAC, 60Hz, 8A
Power Consumption	160W max
Battery Backup	Up to 2hrs