

# ULTRA LOW TEMPERATURE RFID FREEZER TS057

- » REAL-TIME INVENTORY VISIBILITY
- » PROACTIVE TEMPERATURE ALERTING
- » ELIMINATE MANUAL PROCESSES

FORMA  
900 SERIES

# KNOW NOW.

Inventory Management Powered by RFID

## BENEFITS OF THE TS057 - ULT FREEZER

- Maintain temperature range of -85°C to -50°C
- Prevent stock-outs; monitor expiry of products
- Protect critical inventory through secure access
- Provide regulatory compliance
- Tailor with your branding
- Integrate with external temperature monitoring systems



COMPONENT	SPECIFICATION
External Dimensions	78" H x 38" W x 37" D
Internal Dimensions	52" H x 23" W x 19" D
External Space	83" H x 52" W x 81" D Place 6" from back wall Allow 34" for door clearance
Volume	13ft <sup>3</sup>
Weight	725lbs
Cooling Method	2-stage cascade system
Access Control	LF badge

COMMUNICATIONS	SPECIFICATION
RAIN RFID Frequency	902 to 928 MHz
Communications	Ethernet, optional cellular or wifi

POWER	SPECIFICATION
Electrical Power	120VAC, 60Hz, 16A
Power Consumption	1.9KW max
Battery Backup	Up to 2hrs