



TERSO
SOLUTIONS, INC.



LARGE RFID REFRIGERATOR TS074e

- » REAL-TIME INVENTORY VISIBILITY
- » PROACTIVE TEMPERATURE ALERTING
- » ELIMINATE MANUAL PROCESS

knownow

Inventory Management Powered by RFID

BENEFITS OF THE TS074e - LARGE REFRIGERATOR

- Reduce energy consumption by 50% over traditional refrigeration technology
- Maintain temperature range of 2°C to 10°C
- Prevent stock-outs; monitor expiry of products
- Customize product layout
- Protect critical inventory through secure access
- Provide regulatory compliance
- Integrate with external temperature monitoring systems
- Quiet operation provides a more comfortable working environment



| COMPONENT | SPECIFICATION |
|---------------------|-------------------------------------------------------------------------------|
| External Dimensions | 79"H x 30"W x 33"D |
| Internal Dimensions | 58.3"H x 24.8"W x 24.3"D |
| External Space | 87"H x 30"W x 66"D Place 3" from back wall Allow 21" for door clearance |
| Volume | 20.2 ft ³ |
| Weight | 500 lbs. |
| Access Control | Supports LF and HF badge types |
| Warranty | 36 months |

| COMMUNICATIONS | SPECIFICATION |
|------------------|--------------------------------------|
| RFID Frequencies | 902 to 928 MHz |
| Communications | Ethernet, Cellular, Optional WiFi |

| POWER | SPECIFICATION |
|---------------------|-----------------------|
| Electrical power | 115 VAC, 60 Hz, 2.8 A |
| Power Consumption | 3.04 kWh/day |
| RFID Battery Backup | 12 VDC SLA, qty 2 |

www.tersosolutions.com/largerfidrefrigerator

