



UNDERCOUNTER RFID FREEZER TS096

- » REAL-TIME INVENTORY VISIBILITY
- » PROACTIVE TEMPERATURE ALERTING
- » ELIMINATE MANUAL PROCESSES

Inventory Management Powered by RFID

BENEFITS OF THE TS096 - UNDERCOUNTER RFID FREEZER

- Prevent stock-outs; monitor expiry of products
- Transaction-level tracking: know who took what and when
- Customize product layout to best organize your inventory
- Designed to fit under most standard height counter tops
- Protect critical inventory through secure access
- Provide regulatory compliance
- Integrate with external temperature monitoring systems



FEATURE	SPECIFICATION
External dimension (H x W x D)	34.0 in. x 23.75 in. x 24.0 in.
Internal dimension (H x W x D)	24.0 in. x 19.0 in. x 17.5 in.
External clearance (H x W x D)	34.0 in. x 23.75 in. x 25.0 in.
Volume	4.2 ft ³
Unit weight	140.0 lbs.
Drawer capacity	10.0 lbs. per drawer
Operating temperature range	-25.0 °C to -15.0 °C
Temperature uniformity	4 °C
Temperature stability	+/-1.5 °C
Defrost type	Manual
Cooling method	In-shelf tubing
Refrigerant	R600a (SNAP)
Security device	Electronic door latch
Access control	LF/HF Pass Reader

COMMUNICATIONS	SPECIFICATION	ELECTRICAL REQUIREMENTS	SPECIFICATION
RFID frequencies	902 to 928 MHz	Electrical power & frequency	115 VAC, 60 Hz, 2 A
Connection	Ethernet, optional cellular	Power consumption	1.2 kWhr/day
		RFID battery backup	12 VDC SLA, qty 2

*Specifications subject to change without notice
Version 6 9/9/2021