



MIDSIZE RFID FREEZER TS093

- » REAL-TIME INVENTORY VISIBILITY
- » PROACTIVE TEMPERATURE ALERTING
- » ELIMINATE MANUAL PROCESSES

knownow

Inventory Management Powered by RFID

BENEFITS OF THE TS093 - MIDSIZE RFID FREEZER

- Maintain temperature range of -25.0 °C to -15.0 °C
- Prevent stock-outs; monitor expiry of products
- Transaction-level tracking: know who took what and when
- Customize product layout to best organize your inventory
- Protect critical inventory through secure access
- Provide regulatory compliance
- Integrate with external temperature monitoring systems



FEATURE	SPECIFICATION
External dimension (H x W x D)	75.0 in. x 24.0 in. x 24.0 in.
Internal dimension (H x W x D)	63.0 in. x 19.0 in. x 17.0 in.
External clearance (H x W x D)	80.0 in. x 28.0 in. x 47.0 in.
Volume	10.4 ft ³
Device weight	214.0 lbs.
Operating temperature range	-25.0 °C TO -15.0 °C
Temperature uniformity	4 °C
Temperature stability	+/- 1.5 °C
Defrost type	Manual
Cooling method	In-shelf tubing
Refrigerant	R290 (SNAP)
Security device	Electronic door latch
Access control	LF pass reader

COMMUNICATIONS	SPECIFICATION
RFID frequencies	902 to 928 MHz
Connection	Ethernet, optional cellular

ELECTRICAL REQUIREMENTS	SPECIFICATION
Electrical power & frequency	115 VAC, 60 Hz, 3.5 A
Power consumption	2.0 kWh/day
RFID battery backup	12 VDC SLA, qty 2

*Specifications subject to change without notice

