

COMPACT INTERNATIONAL RFID REFRIGERATOR

The Compact International RFID Refrigerator enables automated inventory anywhere in the world with proactive support and 24/7 temperature tracking, plus real-time insight of your inventory. Its compact design includes open shelves or sliding divider trays to suit your inventory organization needs.



FEATURES

AT A GLANCE

Model number	EUR TS060
External dimension H x W x D	1010 mm. x 585 mm. x 620 mm.
Internal dimension H x W x D	618 mm. x 492 mm. x 415 mm.
External Space Required H x W x D	1040 mm. x 685 mm. x 1255 mm. Allow for 50 mm. of clearance on the sides and rear of the device, 30 mm. on the top, and 585 mm. on the front for door clearance
Volume	104 L
Device weight	75 kg
Access control	LF pass reader
Operating temperature range	+2.0 °C to +12.0 °C
Cooling method	Forced air cooled
RAIN RFID frequency	865 to 868 MHz
Communications	Ethernet, optional cellular
Electrical power & frequency	230 VAC, 50 Hz, 2 A
Power consumption	0.7 kWh/day
RFID battery backup	24 V SLA 1.3 Ah (2 count)
Optional riser H x W x D	584 mm. x 559 mm. x 610 mm.
Certifications	RAIN, CE, RoHS



BENEFITS

- Maintain temperature range of +2.0 °C to +12.0 °C
- Prevent stock-outs
- Monitor expiry of products
- Customize product layout
- Protect critical inventory through secure access
- Provide regulatory compliance
- Tailor with your branding
- Integrate with external temperature tracking systems

Compact International
RFID Refrigerator in action

WHAT IS IT GOOD FOR?

APPLICATIONS

**Pharmaceuticals, medications,
and vaccines**
Reagents
Laboratory samples
Refrigerated tissue & biologics
Cardiovascular products
Orthopedic products
**Refrigerated consignment
inventory**

LOCATION USE

Hospitals
Stockrooms
Research labs
Pharmacies
Aerospace
Consignment inventory