

UNDERCOUNTER RFID REFRIGERATOR

Ensure peace of mind with the Undercounter RFID Refrigerator. Leave manual inventory management behind with proactive monitoring and 24/7 temperature tracking, plus real-time insight to your inventory. Its compact design and small footprint allow you to track high-value inventory without compromising space.



FEATURES

AT A GLANCE

Model number	US TS097
External dimension H x W x D	34.0 in. x 24 in. x 24.5 in.
Internal dimension H x W x D	23.5 in. x 17.5 in. x 17.0 in.
External Space Required H x W x D	34.0 in. x 28.0 in. x 52.0 in. Allow for 2 in. of clearance on the sides of the device, 3 in. on the rear, and 24 in. on the front for door clearance
Volume	4.6 ft ³
Device weight	130.0 lbs
Access control	LF pass reader
Cooling method	Forced air cooled
Operating temperature range	+1.0 °C to +10.0 °C
RAIN RFID frequency	902 to 928 MHz
Communications	Ethernet, optional cellular
Electrical power & frequency	120 VAC, 60 Hz, 2 A
Power consumption	1.15 kWh/day
RFID battery backup	24 V SLA 1.3 Ah (2 count)
Optional riser H x W x D	23.0 in. x 22.0 in. x 24.0 in.
Certifications	RAIN, SGS, FCC, SNAP, EnergyStar



BENEFITS

- Maintain temperature range of +1.0 °C to +10.0 °C
- Prevent stock-outs
- Monitor expiry of products
- Customize product layout
- Protect critical inventory through secure access
- Provide regulatory compliance
- Tailor with your branding
- Integrate with third-party inventory temperature monitoring systems

Undercounter RFID Refrigerator in action

WHAT IS IT GOOD FOR?

APPLICATIONS

Pharmaceuticals, medications, and vaccines
Reagents
Laboratory samples
Refrigerated tissue & biologics
Cardiovascular products
Orthopedic products
Refrigerated consignment inventory

LOCATION USE

Hospitals
Stockrooms
Research labs
Pharmacies
Aerospace
Consignment inventory

*Specifications subject to change without notice